

# GC/MS Analysis of Semi-volatile Organic Compounds in Drinking Water: Productivity Solution for US EPA Method 525.2

Jessie Butler, Eric Phillips, Thermo Fisher Scientific, Austin, Texas, USA

## Overview

The analysis of drinking water for semi-volatile organics including thermally labile pesticides of varying polarity is quite challenging. Using a surged splitless injection with the DSQ™ II GC/MS enables laboratories to easily achieve the method detection limits required by United States Environmental Protection Agency (USEPA) Method 525.2. By combining the DSQ II GC/MS with all of the methodologies, consumables, operating procedures, and layered application software, the 525.2 Productivity Solution streamlines method development and validation.

## Results

A mixture of 115 compounds, including polynuclear aromatics, phthalates, Aroclors, and pesticides was analyzed using the DSQ II GC/MS. Over a concentration

range of 0.1 to 10 ng/μL, the calibration curve had an average RSD of 9%. The average method detection limit (MDL) at a 99% confidence level was 0.04 μg/L. Seven replicate injections were made at 5 ng/μL to demonstrate the accuracy of quantitation at a mid-level concentration, which resulted in an average RSD of 1.9% for the compounds, easily meeting the percent recovery of 70 to 130%. US EPA Method 525.2 identifies problem compounds that tend to show higher MDLs, lower recoveries, and potential for non-detection. Table 1 gives an overview of the published MDLs in the EPA Method, the Instrument Detection Limits (IDLs) achieved on the DSQ II, and the calibration curve fit for each problem compound as called out in the method. The EPA Method 525 Productivity Solution produces Gaussian peak shapes and good results for these problematic compounds.

Problem Compounds Listed in Method* with Stated Low Recoveries		MDL Listed in Method (μg/L)	DSQ II Calibration Results		DSQ II IDLs (μg/L)
			% RSD	R <sup>2</sup>	
merphos	breakdown to DEF	Not Detectable (ND)	24	0.9967	0.1320
fenamiphos	poor GC peak	0.95	24	0.9978	0.0320
fenarimol	poor GC peak	1.2	14	0.9993	0.0280
fluridone	poor GC peak	0.55	32	0.9989	0.0490
hexazinone	poor GC peak	0.11	13	0.9991	0.0420
norflurazon	poor GC peak	0.13	19	0.9994	0.0710
stirofos	poor GC peak	0.13	18	0.9997	0.0140
tebuthiuron	poor GC peak	2.8	16	0.9993	0.0300
tridemefon	poor GC peak	0.33	11	0.9998	0.0420
tricyclazole	poor GC peak	2.6	26	0.9959	0.0610
anthracene	photodegradation	0.068	3	0.9995	0.0140
benzo(a)anthracene	photodegradation	0.2	3	0.9993	0.0270
dimethyl phthalate	in lab at ~ 2 μg/L	0.058	3	0.9999	0.0270
atraton	ND at pH of 2	0.16	10	0.9998	0.0240
prometon	ND at pH of 2	0.38	10	0.9999	0.0540
carboxin	ND in water	1.4	16	0.9995	0.0270
disulfoton	ND in water	0.62	8	0.9999	0.0270
metribuzin	< 50%	0.16	8	0.9999	0.0200
cyanazine	ND at pH of 2	0.17	10	0.9999	0.0580
pentachlorophenol	loss of sensitivity < 5 ng	ND	14	0.9948	0.4900

Table 1: EPA Method Problem Compounds, published MDLs, and DSQ II Instrument Detection Limits (IDLs)

\*Characteristics as described in the method in Section 13.2 Problem Compounds<sup>1</sup>

### Keywords:

- DSQ II GC/MS
- Drinking Water
- EnviroLab Forms 2.0 Software
- Semi-Volatile Analysis
- Splitless Injection Method
- US EPA Method 525.2

## Methods

In full-scan mode, the DSQ II offers excellent sensitivity, which allows a conventional splitless injection to be used instead of more difficult techniques, such as cold on-column or temperature programmed injections. A 1  $\mu$ L injection was made in the splitless mode with a surge pressure of 37 psi and a 1 mL/min column flow. A proprietary phase, 30 m x 0.25 mm i.d. x 0.25  $\mu$ m film thickness TRACE™ TR-525 column was used (Thermo Scientific PN 26RX142P). The column and inlet were tested for activity by injecting the performance mix, a solution of DFTPP, endrin, and 4,4'-DDT at 5  $\mu$ g/mL in iso-octane (Thermo Scientific PN 60181-311). The mass spectrometer was tuned automatically using Autotune Tuning software to meet DFTPP tuning criteria. The TIC for the mid-level standard is shown in Figure 1.




Figure 1: TIC of mid-level standard (5 ng/ $\mu$ L)

Data was acquired using Xcalibur™ instrument control software and then processed automatically using the EnviroLab™ Forms 2.0 software. Complete methods are available for EPA Methods and reports can be saved in \*.pdf, \*.doc, \*.xls, or \*.rtf format. Reports can also be uploaded to a laboratory information management system (LIMS) using \*.XML format files.

## Conclusion

This method provides fast chromatography while maintaining optimal separations of the chromatographic peaks. The DSQ II successfully met the QC criteria for EPA Method 525.2 in a splitless mode using the TRACE GC Ultra™. Combined with a fast scanning rate, this Productivity Solution features excellent separation and sensitivity, enabling the generation of MDLs required in the EPA method. The Productivity Solution (Figure 2) contains all of the consumables needed to perform the instrument validation, including standard solutions and vials. The electronic methods, Standard Operating Procedure, and How To Manual with Quick Start Guide can be automatically downloaded to the local computer by installing the Interactive Reference CD (PN 120296-CD). A Validation Data CD (PN 120296-VCD) may be used as a reference of typical data for the method.



Figure 2: EPA Method 525.2 Semi-volatiles in Water Productivity Solution

EnviroLab Forms 2.0 software matches the workflow in environmental laboratories around the world, is simple to use, and allows novice users to be instantly productive. Now, adding a DSQ II system to an existing laboratory workflow is easier than ever.

For detailed information on the instrument and processing parameters, as well as calculated values for all compounds analyzed, please visit our website at [www.thermo.com/gc](http://www.thermo.com/gc) and request TN10172.

## References

1. Method 525.2 Determination of Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column Gas Chromatography/Mass Spectrometry, Rev 2.0, National Exposure Research Laboratory Office of Research and Development U. S. Environmental Protection Agency Cincinnati, Ohio 45268

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

**Australia**  
+61 2 8844 9500  
**Austria**  
+43 1 333 50340  
**Belgium**  
+32 2 482 30 30  
**Canada**  
+1 800 532 4752  
**China**  
+86 10 5850 3588  
**Denmark**  
+45 70 23 62 60  
**France**  
+33 1 60 92 48 00  
**Germany**  
+49 6103 408 1014  
**India**  
+91 22 6742 9434  
**Italy**  
+39 02 950 591  
**Japan**  
+81 45 453 9100  
**Latin America**  
+1 608 276 5659  
**Netherlands**  
+31 76 587 98 88  
**South Africa**  
+27 11 570 1840  
**Spain**  
+34 91 657 4930  
**Sweden/Norway/Finland**  
+46 8 556 468 00  
**Switzerland**  
+41 61 48784 00  
**UK**  
+44 1442 233555  
**USA**  
+1 800 532 4752

[www.thermo.com](http://www.thermo.com)



Thermo Electron Corporation,  
Austin, TX USA is ISO Certified.

©2006 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

AN10173\_E 12/06M